

# OncoStrands™

## Extended Panel (Tissue Biopsy)

OncoStrands™ Extended Panel (Tissue) tests for 108 selected cancer-relevant genes in a pan-cancer setting from formalin-fixed, paraffin embedded (FFPE) or cytology tumour material. The test was designed with poor quality FFPE samples in mind, and it provides high accuracy, sensitivity, and specificity for even heavily degraded FFPE tissue samples.

In contrast to hotspot panels, this hybrid capture next generation sequencing assay covers the full coding region of 103 out of 108 genes, screening for mutations, small indels, promoter mutations (*TERT* gene), and copy number variations (CNVs). The assay is particularly useful for advanced pancreatic, prostate, breast, hepatobiliary, colorectal, uterine, ovarian and fallopian tube tumours as most of the currently known DNA damage repair genes, including *BRCA1* and *BRCA2* are covered. It also provides an accurate microsatellite instability (MSI) value for each sample (close to 170 MSI loci have been included in the panel design compared to just 5 loci screened for in conventional tests).

The assay can detect biomarkers that are frequently mutated in various cancer types. This data, combined with cutting-edge curation and analysis solutions, enables not only the detection of biomarkers but also matching with key and emerging treatment guidelines, FDA-approved targeted therapies, and clinical trials. The assay provides an accurate MSI score, and even more value when combined with the 18 gene OncoStrands™ Essential RNA fusion panel.

This enables oncologists to select the most appropriate therapeutic approach, anticipate prognosis of the disease course, and fully personalize the disease management for each patient.



**108 Pan-Cancer Gene DNA Targets**  
SNVs, CNVs, MSI  
Full coding regions – 103 genes including *BRCA1/2*



**Low (20%)\* Tumour Content (TC) Required**



**Screens Key and Emerging Biomarkers for Multiple Tumour Types**  
Screens for FDA-approved therapies & NCCN recommended biomarkers



**High Accuracy at ≥5% Limit of Detection**  
Customised assay chemistry validated\* even for poor quality FFPE samples



**10 Days Turnaround Time\*\***



**Bespoke Consultation with Molecular Pathologists**

\*Validated at 10% TC for SNVs and 20% TC for CNVs and MSI

\*\*Comprehensive report available 10 working days from laboratory sample receipt and subject to sample acceptance criteria

### Test Specifications & Validation Characteristics

Based on in-house validation of clinical samples and reference standards

<b>Methodology</b>	Next generation sequencing
<b>Aberrations Covered</b>	SNVs, indels, CNVs, MSI
<b>Specimen Requirements</b>	<ul style="list-style-type: none"> <li>FFPE tissue block OR minimum 10 unstained sections (each 5µm thick). Minimum tumour content of 20%.</li> <li>Copy of histology report.</li> </ul>

\*Please refer to the Molecular Oncology Request Form for full specimen requirements.

Mutation Type	Accuracy	Sensitivity	Specificity	Limit of Detection
SNVs/indels	100%	100%	100%	≥5%
CNVs*	100%	99.0%	100%	N/A
MSI	MSI -percentage concordance with samples run on orthogonal tests = 100% (close to 170 MSI loci included for screening in the assay)			

\*CNVs on targeted NGS platforms is an estimate based on prediction algorithm which considers multiple factors. The assay is validated for gene amplifications of ≥5 and homozygous deletions.

### OncoStrands™ Extended Panel Full Gene List

*AKT1, ALK, APC, AR, ARID1A, ATM, ATR, BAP1, BARD1, BRAF, BRCA1, BRCA2, BRIP1, CDH1, CD274, CDK4, CDK6, CDK12, CDKN2A, CDKN2B, CHEK1, CHEK2, CSF1R, CTNNA1, CYSLTR2, DDR2, EGFR, EPCAM, ERBB2, ERBB3, ERCC2, ESR1, EZH2, FANCA, FANCL, FBXW7, FGFR1, FGFR2, FGFR3, FOXL2, GNA11, GNAQ, GNAS, HDAC2, HNF1A, HRAS, H3F3A(hotspots only), H3F3B(hotspots only), HIST1H3B(hotspots only), HIST1H3C(hotspots only), IDH1, IDH2, JAK2, KEAP1, KIT, KRAS, MAP2K1, MAP2K2, MET, MLH1, MRE11, MSH2, MSH6, MYC, MYCN, MTOR, MUTYH, NBN, NF1, NF2, NFE2L2, NOTCH1, NRAS, NTRK1, NTRK2, NTRK3, PALB2, PDGFRA, PDGFRB, PIK3CA, PMS2, POLD1, POLE, PPP2R1A, PPP2R2A, PTCH1, PTEN, PTPN11, RAD50, RAD51B, RAD51C, RAD51D, RAD54L, RAF1, RB1, RET, ROS1, SETD2, SF3B1, SMAD4, SMARCB1, SMO, SRC, STK11, TERT(5' promoter only), TP53, TSC1, TSC2, and VHL*



## This Assay Screens For Genetic Alterations That Are Linked To Current FDA Approved Therapies\*

Tumour Type	Genetic Alterations/Biomarker	FDA-Approved Targeted Therapies (OR Contraindications)
<b>FDA-Approved Treatments For Specific Genetic Alterations In Specific Tumour Types</b>		
Non-small-cell lung cancer (NSCLC)	BRAF V600E	Dabrafenib + Trametinib
Melanoma		Dabrafenib, Vemurafenib Dabrafenib + Trametinib, Encorafenib + Binimetinib, Vemurafenib + Cobimetinib, Trametinib
Anaplastic thyroid cancer		Dabrafenib + Trametinib
Colorectal cancer		Encorafenib + Cetuximab
Melanoma	BRAF V600K	Dabrafenib + Trametinib, Encorafenib + Binimetinib, Vemurafenib + Cobimetinib, Trametinib
Ovarian cancer, fallopian tube cancer, peritoneal cancer	Deleterious or suspected deleterious germline or somatic mutations in BRCA1 and/or BRCA2	Olaparib, Rucaparib, Niraparib
Prostate cancer		Olaparib, Rucaparib
Ovarian cancer, pancreatic adenocarcinoma	Deleterious or suspected deleterious germline mutations in BRCA1 and/or BRCA2	Olaparib
HER-2 negative breast cancer		Olaparib, Talazoparib
Prostate cancer	Deleterious or suspected deleterious germline or somatic mutations in ATM, BARD1, BRIP1, CDK12, CHEK1, CHEK2, FANCL, PALB2, RAD51B, D51C, RAD51D, and RAD54L	Olaparib
NSCLC	EGFR exon 19 deletions, L858R	Afatinib, Dacomitinib, Erlotinib, Gefitinib, Osimertinib
	EGFR exon 20 insertions	Amivantamab
	EGFR non-resistant mutations other than exon 19 deletions and L858R	Afatinib
	EGFR T790M	Osimertinib
Breast cancer	ERBB2 amplification	Ado-Trastuzumab Emtansine, Capecitabine + Trastuzumab + Tucatinib, Neratinib, Pertuzumab + Trastuzumab, Trastuzumab, Trastuzumab Deruxtecan
		Trastuzumab
		Trastuzumab Deruxtecan
Oesophageal cancer		
Gastric cancer, gastroesophageal junction cancer		
Bladder cancer	FGFR3 oncogenic mutations	Erdafitinib
NSCLC	KRAS G12C	Sotarasib
NSCLC	MET exon 14 skipping	Capmatinib, Tepotinib
Colorectal cancer	MSI-H	Ipilimumab + Nivolumab, Nivolumab
Endometrial cancer		Dostarlimab
Gastrointestinal stromal tumour (GIST)	PDGFRA exon 18 mutations	Avapritinib
HR+ HER2- breast cancer	Oncogenic mutations in PIK3CA	Fulvestrat + Alpelisib
Medullary thyroid cancer	Oncogenic mutations in RET	Pralsetinib, Selpercatinib
<b>FDA-Approved Treatments For Specific Biomarkers In Tumour Type-agnostic Indications</b>		
Pan-cancer (solid tumours)	MSI-H	Pembrolizumab
<b>FDA-Approved Treatments That Are Not Biomarker-linked In Solid Tumours Characterised By Specified Genetic Alterations</b>		
Neurofibroma	Oncogenic mutations in NF1	Selumetinib
Epithelioid sarcoma	SMARCB1 deletions	Tazemetostat
Subependymal giant cell astrocytoma (SEGA)	Oncogenic mutations in TSC1/TSC2	Everolimus
GIST	KIT exon 9, 11, 13, 14, 17 mutations	Imatinib, Sunitinib (post progression on Imatinib), Regorafenib (post progression on Imatinib and Sunitinib), Ripretinib (post progression on ≥3 kinase inhibitors including Imatinib)
<b>FDA-Listed Genetic Alterations Contraindicated For Specific Treatments</b>		
Colorectal cancer	KRAS and/or NRAS exon 2, 3, and 4 mutations	Contraindicated for Panitumumab, Cetuximab
Pan-cancer (solid tumours)	NTRK1 and NTRK3 known acquired resistance mutations (e.g., NTRK1 G595R and G667C; NTRK3 F617L, G623R, and G696A)	Contraindicated for Entrectinib, Larotrectinib
<b>FDA-Approved Combination Treatments With Nontargeted Therapies For Specific Genetic Alterations</b>		
Melanoma	BRAF V600	Atezolizumab + Cobimetinib + Vemurafenib
Fallopian tube, ovarian, primary peritoneal carcinoma	Deleterious germline or somatic mutations in BRCA1 and/or BRCA2	Bevacizumab + Olaparib
NSCLC	EGFR exon 19 deletions, L858R	Erlotinib + Ramucirumab
Esophagogastric cancer	ERBB2 amplification	Trastuzumab + Cisplatin + Capecitabine or Fluorouracil

\*Modified from www.ascopubs.org, current as of June 2022

### Reports Include

- Contents as per the latest AMP and CAP guidelines
- Recommended clinical matching with biomarkers, and clinical trials as per FDA, EMA, NCCN, ESMO, etc.

### Services Include

- Quality control for tissue adequacy performed by staff pathologist
- Tests run in house by qualified scientific and clinical staff under an accredited environment
- Complimentary consultation on various aspects of testing (e.g., appropriate test options based on tumour type, tissue availability etc.) provided by qualified staff molecular pathologist

### Additional Services

- IHC- MMR, PDL-1, ALK, ROS1
- Range of Oncostrands™ (oncosomatic) and hereditary panels

### About Us

At LifeStrands Genomics laboratories we believe that everyone should have access to better healthcare through the advancement of clinical genomics. Within our accredited laboratories, our dedicated team of medical professionals and scientists work together to deliver high-quality and reliable genomic solutions to clinicians, patients & researchers.



## LifeStrands Genomics Laboratory Locations

### Australia

Building A (G.01) 18-24 Ricketts Road, Mount Waverley, Victoria 3149, Australia

W www.lifestrandsgx.com.au  
E genomics@lifestrandsgx.com

### Singapore

The Galen 61 Science Park Road #03-13/14, Singapore 117525

W www.lifestrandsgx.com  
E enquiry@lifestrandsgx.com

### Malaysia

No. 47-3A, Level 3, Jalan PJU 5/12 Dataran Sunway, Kota Damansara 47810 Petaling Jaya, Selangor, Malaysia

W www.lifestrandsgx.com  
E enquiry@lifestrandsgx.com